

COUNTLESS
POSITION
COMBINATIONS



Shown In Ultraleather™
Fudge Fabric

INFINITY COLLECTION

Enjoy endless relaxation with the **Infinity Collection** from **VivaLift!™** Power Lift Recliners. True-Infinite 4-motor recliner with power lumbar offers personalized comfort, and the full width power headrest lifts your head, neck and shoulders. The Infinity Collection offers true infinite positioning and a footrest extension that lets you truly stretch out and relax. A convenient hand control with USB charger port lets you charge your smart device, and the industry-first Lithium battery backup ensures your lift chair functions when the power goes out. With an array of gorgeous fabrics, it's an easy decision to make this fabulous collection a part of your home.

**THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!**

(US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com

Pride®
Mobility Products Corp.
— LIVE YOUR BEST! —

INDEPENDENT FUNCTIONS:



 TRENDELENBURG



 POWER HEADREST



 POWER LUMBAR



 BACKREST




 FOOTREST EXTENSION
(Standard)



 UPWARD LIFT

STANDARD FABRIC SELECTION GUIDE

DURASOFT 



OAT



TIMBER

OPTIONAL INFINITY FABRIC SELECTIONS

ULTRALEATHER™  



BUFF






FUDGE

SPECIFICATIONS

Model	PLR-525iM
Position	True Infinite-Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 6'
Weight Capacity	375 lbs.
Seat-to-Floor	20"
Seat Depth	20"
Seat Width	20"
Top of Back to Seat	26"
Back Style	Blown fiber, Split back
Fully Padded Chaise	Yes
Footrest Extension	Standard
Distance From the Wall	20"
Overall Width	33"
Chair Weight	148 lbs.
Warranty	7-Yr. Pro-rated

STANDARD FEATURES

-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  TRUE-INFINITE POSITION

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

(US) 800-800-8586 | (CAN) 888-570-1113
 pridemobility.com | prideliftchairs.com

Pride
 Mobility Products Corp.
 — LIVE YOUR BEST —