

An Active Lifestyle Deserves a Sporty Ride



Shown in
Candy Apple Red

VICTORY[®] SPORT



Shown in
Viper Blue

The **Victory[®] Sport** incorporates a sporty design with extra comfort and performance features for a truly unique ride experience. With a maximum speed of 8 mph—the fastest within the Victory Series—this scooter will accelerate your active lifestyle. Standard features include feather touch disassembly, wraparound delta tiller, 10" solid tires, front and rear suspension, high-back reclining seat, directional signals, a high-intensity LED lighting package, an easy-to-read LED battery meter, and a front basket. Available in Viper Blue and Candy Apple Red.

PRIDE[®] MOBILITY SCOOTERS
LIVE YOUR BEST[®]

pridemobility.com





Specifications 4W:SC710 DXW

Weight Capacity	350 lbs.
Maximum Speed ¹	Up to 8 mph
Ground Clearance ³	2.25" at motor, 1" at anti-tip, 3" at rear lock up
Turning Radius ³	51.75"
Overall Length ³	47"
Overall Width ³	26"
Seat-to-Ground Height	21.5"-24.25"
Seat-to-Deck Height	16"-18.5"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	Front and rear
Range Per Charge ^{1,2}	Up to 11 miles
Total Weight Without Batteries ^{3,7}	173 lbs.
Heaviest Piece When Disassembled ³	68 lbs. (front section)
Standard Seating	Type: Reclining high-back with sliders Weight: 45.5 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Electronic, regenerative and electromechanical
Battery Requirements ^{3,4,5}	Size: (2) 40 AH Weight ⁵ :32 lbs. each
Battery Charger	3.5-amp, off-board (U-1), 5-amp, off-board (40 AH)
Colors	 

Candy Apple Red Viper Blue

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- AGM or gel cell type required.
- Battery weight may vary based on battery manufacturer.
- Includes standard seat weight.

Features

- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Delta tiller with wraparound handles
- Exclusive stylish, lightweight, non-scuffing, black, low-profile tires
- Easy-to-read LED battery meter
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Black vinyl reclining high-back seat with sliders, height adjustable head rest, flip-up, height and width adjustable armrests
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port

Accessories

- Cane/crutch holder (3" single, 3" double)
- Cell phone holder
- Cup holder
- Double crutch holder
- Forearm crutch holder
- Lap belt (50", 60", 70")
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

Optional Seating

- Optional 22" x 18" high-back seat

Military Patches



LIFT FINDER

Pride offers a variety of lifts to **transport your mobility product with ease**. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

(US) 800-800-8586

(CAN) 888-570-1113

pridemobility.com

