


Bobcat X4 Transportable Scooter



INSET A




INSET B



INSET C



INSET D



INSET E



INSET F

- Ultra-lightweight for superior portability
- Quick and easy disassembly
- Color-through shroud hides scuffs and scratches
- Angle adjustable flip-up armrests for maximum comfort and convenient boarding and exiting
- Easy finger-pull throttle control
- Height adjustable seat with fold down backrest
- Flat-free, non-marking tires for worry-free travel
- Large, metal carry basket
- Easy-to-adjust tiller angle places the scooter's controls in the ideal position
- Anti-tip wheels for additional safety
- Easy freewheel operation

PRODUCTS

PRODUCTS	Description
BOBCATX4	Red, BOBCAT X4 4-WHEEL SCOOTER RED, 1/ea

SPECIFICATIONS

SPECIFICATIONS	Description
Weight Capacity	265 lbs.
Top Speed	4 mph
Maximum Range	7.5 miles
Turning Radius	53"
Climbing Angle	6°
Ground Clearance	2.5"
Dimensions (L x W)	42" x 19.25"
Floor-to-Seat Height	17"-19"
Seat Dimensions (W x D)	16.5" x 13.5"
Controller	PG S-Drive 45A
Motor	24V x 270W x 4700 rpm
Batteries (pair)	12V x 10AH
Charger	1.5A Offboard
Brakes	Electromagnetic
Wheels	8" x 2" Flat-Free
Base Weight	59 lbs.
Battery Weight (pair)	18 lbs.
Seat Weight	16 lbs.
Total Weight	93 lbs.
Heaviest Piece Weight (Front Section)	38 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components	24 Months

SPECIFICATIONS

Description

Warranty on Batteries

12 months

ACCESSORIES

Description

SF8010

Oxygen Tank Holder, 1/ea

SF8020

Rear Basket, 1/ea

SF8030

Cane/Crutch Holder, 1/ea

SF8050

Sun Shade, 1/ea

PARTS

Description

BC31009N

C12-071-00100

BC31208B

C08-050-00301

C18-071-00100

C18-071-00200

C08-071-00100

C02-071-00100

S31003

BC35051

C08-069-01001

C02-069-00401

C29-016-00301

C01-069-00103

C01-071-00100

C08-069-04800

S03-069-00200

C03-069-00101

BC350103

C03-010-00305

BC31142

BC35046

S31026

C08-071-00600

C08-071-00200

C08-071-00400

C08-071-00512

C01-070-00300

BC350112

C09-071-00200

C09-071-00300

BC31132

S350016910

C18-071-00601

C18-071-00800

BATT12V10AH

S07-071-00100

C07-071-00400

C07-071-00300

C07-071-00301

S31052

C07-071-00100

C07-071-00101

C18-071-00300

BC40AFUSE

C07-050-00203-C