

Spitfire EX2 4-Wheel Travel Scooter with 12AH Battery



WHITE



BLUE



SILVER



BLACK



PINK



INSET F



INSET G



INSET H

- 6 sets of interchangeable color panels in red, blue, silver, black, white and pink
- Swivel mount cup holder included
- Extra bright LED headlight
- Charge ports located on the tiller and battery pack for easy charging
- Height adjustable swivel seat
- Adjustable padded armrests with reflectors
- Large plastic basket with removable liner
- Disassembles into 5 easy to carry pieces
- Delta tiller for easy throttle control

PRODUCTS	Description
SFEX2417FS-12	KO800, 1/ea
SFEX2417FS-21	KO800, 1/ea

SPECIFICATIONS	Description
Weight Capacity	300 lbs.
Max Speed	5 mph
Maximum Range	9 Miles 15 Miles
Max. Climbing Angle	6°
Turning Radius	54"
Ground Clearance	3"
Overall Dimensions	42" x 20" x 36"
Seat Width	17"
Seat Depth	16"
Back Height	16"
Seat-to-Floor Height	19" – 21"
Motor	24V x 270W x 4700 rpm
Controller	45A
Batteries	12AH 21AH
Battery Charger	2A Offboard
Brakes	Electromagnetic
Freewheel Mode	Yes
Front Wheels	7.5" x 2"
Rear Wheels	8.5" x 3"
Base Weight	66 lbs.
Battery Weight	22 lbs. (12AH) 28 lbs. (21AH)
Seat Weight	17 lbs.
Total Weight	105 lbs. (12AH) 111 lbs. (21AH)

SPECIFICATIONS

Description

Warranty on Frame	Limited Lifetime
Warranty on Electronics	24 Months
Warranty on Batteries	12 Months

ACCESSORIES

Description

SF8010	Oxygen Tank Holder, 1/ea
SF8020	Rear Basket, 1/ea
SF8030	Cane/Crutch Holder, 1/ea