


Panther Heavy Duty 4-Wheel Scooter



INSET A



INSET B



INSET C



INSET D

- Front and rear suspension on 14" pneumatic tires for the smoothest ride possible
- One-touch, infinite tiller adjustment allows the user to easily position the tiller into the ideal driving position
- Extra-bright LED lighting for superb visibility
- Ergonomic delta control handles for ease of steering and throttle control
- Large front plastic basket with carry handle
- 20" – 22" wide Captain's Seat is height adjustable and features a fold-down, reclining backrest with adjustable headrest
- Hazard lights, turn signals and reflectors for superb visibility
- Dual rear view mirrors
- Positioning belt included
- Anti-tip wheels
- Easy, freewheel mode
- USB port for charging your mobile devices

PRODUCTS

Description

PANTHER20CS	Panther HD Scooter 20 Capt.Seat, 1/ea
PANTHER22CS	Panther HD Scooter 22 Capt.Seat, 1/ea

SPECIFICATIONS

Description

Weight Capacity	425 lbs.
Top Speed	8 mph
Estimated Range	25 miles
Max Climbing Angle	7.5°
Turning Radius	53"
Ground Clearance	4"
Overall Dimensions	56" (L) x 26" (W) x 53" (H)
Seat Dimensions	20-22" (W) x 20" (D) x 22" (H)
Motor	800W x 4500RPM
Controller	Dynamic Rhino 2 120A
Batteries	50AH
Battery Charger	5A
Brakes	Electromagnetic
Freewheel Mode	Yes
Front Wheels	14" x 4"
Rear Wheels	14" x 4"
Anti-Tip Wheels	2.5"
Base Weight	145 lbs.
Battery Weight (Pair)	64 lbs.
Seat Weight	55 lbs.

SPECIFICATIONS	Description
Total Weight	264 lbs.
Warranty	24 Months on Electronics 12 Months on Batteries
Warranty on Frame	Limited Lifetime