


ZooMe 3-Wheel Recreational Scooter



INSET A



INSET B



INSET C



INSET D



INSET E

- Powerful hub driven motor providing speeds up to 15 mph
- Folding tiller for easy transportation in a vehicle, or for storage
- 2 keys included
- Hand brake with locking mechanism for parking on inclines
- Forward and reverse controls are easily changed with the switch of a button
- Fast and slow speed settings easily changed with the switch of a button
- Front basket for storage
- Pneumatic tires for the smoothest ride possible
- Extra bright LED headlight
- Anti-tip wheel

PRODUCTS

Description

ZOOME3

Red, 1/ea

SPECIFICATIONS

Description

Weight Capacity	300 lbs.
Top Speed	15 mph
Estimated Range	20 miles
Maximum Climbing Angle	10°
Turning Radius	22"
Ground Clearance	4.25"
Overall Dimensions	41" (L) x 27" (W) x 28.5" (H)
Seat Dimensions	14.5" (W) x 10" (D)
Motor	350W Brushless Hub
Controller	350W
Batteries	3ea 12V x 12AH
Battery Charger	Offboard Smart Charger
Brakes	Hand Brakes
Freewheel Mode	No
Front Wheels	16" x 2"
Rear Wheels	4" x 2"
Anti-Tip Wheels	Standard
Base Weight	78 lbs.
Warranty	3 Years on Main Frame, 12 Month on Electrical Harness, Controller, Motor/Gearbox, 6 Months on Batteries