

LiteRides 3 wheel GL111D 4 wheel GL141D

Seat Color Black Vinyl

Body Colors Red BI

Enjoy the best of both worlds with the LiteRider™ scooter! The LiteRider is a full-size scooter with all of the benefits of a compact, portable scooter! The LiteRider comes standard with a Delta Tiller, a large 17"x16" stadium style seat, and adjustable armrests. With a 300 lb. weight capacity, the LiteRider can be quickly and wirelessly disassembled for transport or storage.





Model	GL111D	CLAMAD
Medicare Code		GL141D
	Not Coded 300 lbs	Not Coded 300 lbs
Weight Capacity Drive Wheels		Rear
	Rear 5 mph	
Maximum Speed Operating Range up to 1	12 miles	5 mph 12 miles
Ground Clearance to center deck	3.75"	3.75"
Ground Clearance to transaxle/		
motor	2.5"	2.5"
Turning Radius	39"	48"
Type Batteries	2-22AH	2-22AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
3	90 Amp	90 Amp
Controller	Dynamic R	Dynamic R
	Series	Series
Dimensions		
Length	41.5"	43.25"
Width	22"	22"
Height (ground to top of back of seat)	37" - 39"	37" - 39"
Ground to Top of Seat	23.5" - 25.5""	23.5" - 25.5"
Ground to Top of Deck	6"	6"
Deck to Top of Seat	18" - 20"	18" - 20"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	30"	30"
Weight of Unit: (assembled)	123 lbs	131 lbs
Front Half	44 lbs	52 lbs
Rear Half	28 lbs	28 lbs
Seat w/Arms	22 lbs	22 lbs
Batteries	29 lbs/pack	29 lbs/pack
Tires:	Solid/Low Profile	Solid/Low Profile
Front	9"	9"
Rear	9"	9"
Rear Anti-Tippers	2"	2"
Standard Seat	Stadium Style	Stadium Style
Back Height - no headrest	16"	16"
Width x Depth	17" X 16"	17" X 16"
Color	Black Vinyl	Black Vinyl
Warranty:		,
Frame	Lifetime	Lifetime
Drive Train	2 Years	2 Years
Electronics <sup>2</sup>	2 Years	2 Years
Batteries	1 Year	1 Year
<del></del>		541

- 1. Battery range tested at 200 lbs., and will vary due to rider weight, drive surface and terrain.
- <sup>2.</sup> Electronics warranty excludes batteries.

## **Features**

- **Delta Tiller standard**
- 300 lb. weight capacity
- Rugged anti-rust steel frame tubing
- Rugged plastic floor deck for optimum durability
- Disassembles quickly and easily

## Enjoy quick and wireless disassembly!







Pull out the battery pack



Remove the baskets from under the seat and from the tiller (not shown)



To separate the rear transaxle from the frame, squeeze the release lever located on the underside of the frame



While squeezing the release lever, pull up on the rear frame



Disassembles into six pieces for transport and storage! No wires to connect or disconnect!



Brand new LED ultra bright adjustable headlight





