







## **Escape**

The Escape scooter is designed for the user on the go. It has a 250 lb. weight capacity and has a quick-disconnect system for effortless disassembly. The unit disassembles into three parts, the heaviest weighing just 33 lb. The Escape maneuvers with ease and has an adjustable tiller for added comfort during operation. The Escape also comes standard with solid tires for a smoother ride and swivel seating system. The Escape also features a split battery pack design that make it convenient for travel.

- + 250lb Weight Capacity
- + Heaviest Part Weighs 33lb
- + Dismantles into Three Parts
- + Swivel Seating
- + Curved Tiller



Red



Blue

Model Number: 7A

## **Technical Specifications**

| Overall Dimensions (L x W x H)       | in                     | 43.7" x 21.3" x 36.4" |
|--------------------------------------|------------------------|-----------------------|
| Number, Size of Tires                | Front, Rear            | 2, 8" (2, 8")         |
| Suggested User Weight <sup>1</sup>   | lb                     | Max. 250              |
| Weight of Heaviest Part <sup>2</sup> | lb                     | 33                    |
| Weight <sup>3</sup>                  | lb                     | 110                   |
| Each Battery Module Weight           | lb                     | 10                    |
| Battery Size                         | Type x PCS             | 12V14AH x 2           |
| Charger                              | Amp                    | Off board (2A)        |
| Output Power of Motor                | HP                     | 0.45                  |
| Maximum Speed                        | MPH                    | 4                     |
| Maximum Climbing Angle/              | %(°)                   | 10(6)                 |
| Safe Climbing Angle <sup>4</sup>     |                        |                       |
| Range (Per Charge) <sup>6</sup>      | Mile                   | 8.8                   |
| Turning Radius                       | in                     | 46"                   |
| Ground Clearance                     | in                     | 1.5"                  |
| Handle Bar                           | Туре                   | T-Bar                 |
| Model #                              |                        | 7A                    |
| Standard Seat                        | $W \times D \times BH$ | SF181614B             |

**Note:** <sup>1</sup> On level Road. <sup>2</sup> Without basket and seating. <sup>3</sup> With basket, seat & batteries. <sup>4</sup> Base on user weight of 200LB. <sup>5</sup> Based on user weight of 165LB. <sup>6</sup> After the batteries and mechanical moving parts fully break in.