



Scootie

The Scootie is light-weight, easy to transport and has articulating front wheels. The Scootie weighs in at just 82 lb. with the batteries installed. The patented articulating front wheels allow four-wheel stability with three-wheel maneuverability. The Scootie is length adjustable and can fold down for easy transport and storage.

- + 250lb weight capacity
- + Length Adjustable
- +Removable Battery Pack
- + Keyed Ignition Security
- + Adjustable Flip-up Arms



Blue

Model Number: TE-787NA Medicare Billing Code: K0800

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

34"

Technical Specifications

Overall Dimensions (L x W x H)	in	37-40″ x 17″ x 34
Number, Size of Tires	Front, Rear	2, 7" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	25
Weight ³	lb	82
Each Battery Module Weight	lb	24
Battery Size	Type x PCS	12V14AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.75
Maximum Speed	MPH	4
Maximum Climbing Angle/	%(°)	10(6)
Safe Climbing Angle ⁴		
Range (Per Charge) ⁶	Mile	10
Turning Radius	in	32″
Ground Clearance	in	1.5″
Handle Bar	Туре	T-Bar
Model #		TE-787NA
Standard Seat	W x D x BH	SF151413BG

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.

