







Sunrunner 3

The Sunrunner 3 is a mid-size three-wheel scooter with a 300 lb. weight capacity. The unit comes standard with a deluxe captain's seat that can accommodate most users. It has an adjustable tiller for added user comfort and a front head light for added convenience for evening operation. It is available optionally with a full lighting package adding turn indicators and hazard lights. The unit has a 3in. ground clearance for operation on multiple terrains.

- + 300lb Weight Capacity
- + Slide & Swivel Seating
- + Adjustable Tiller
- + Optional Rear Light Package
- + Adjustable Seat Height



Burgundy



Blue

Model Number: 888B-3 Medicare Billing Code: K0806

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

| Overall Dimensions (L x W x H) | in | 44" x 23" x 42" |
|--------------------------------------|------------------------|-----------------|
| Number, Size of Tires | Front, Rear | 1, 9" (2, 10") |
| Suggested User Weight ¹ | lb | Max. 300 |
| Weight of Heaviest Part ² | lb | _ |
| Weight ³ | lb | 165 |
| Each Battery Module Weight | lb | 24 |
| Battery Size | Type x PCS | 12V35AH x 2 |
| Charger | Amp | Off board (5A) |
| Output Power of Motor | HP | 1 |
| Maximum Speed | MPH | 5 |
| Maximum Climbing Angle/ | %(°) | 14(8) |
| Safe Climbing Angle ⁴ | | |
| Range (Per Charge) ⁶ | Mile | 25 |
| Turning Radius | in | 37" |
| Ground Clearance | in | 3" |
| Handle Bar | Туре | T-Bar |
| Model # | | 888B-3 |
| Standard Seat | $W \times D \times BH$ | SV191821B |

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.