




## Sunrunner 4

The Sunrunner 4 is a mid-size four-wheel scooter with a 300 lb. weight capacity. The unit comes standard with a deluxe captain's seat that can accommodate most users. It has an adjustable tiller for added user comfort and a front head light for added convenience for evening operation. It is available optionally with a full lighting package adding turn indicators and hazard lights. The unit has a 3in. ground clearance for operation on multiple terrains.

- + 300lb Weight Capacity
- + Slide & Swivel Seating
- + Adjustable Tiller
- + Optional Rear Light Package
- + Adjustable Seat Height



Model Number: 888B-4 Medicare Billing Code: K0899

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

## **Technical Specifications**

48" x 23" x 43" Overall Dimensions (L x W x H) in Front, Rear Number, Size of Tires 2, 10" (2, 10") Suggested User Weight<sup>1</sup> lb Max. 300 Weight of Heaviest Part<sup>2</sup> lb Weight <sup>3</sup> lb 180 **Each Battery Module Weight** lb 24 **Battery Size** Type x PCS 12V35AH x 2 Off board (5A) Charger Amp ΗP **Output Power of Motor** 1 **Maximum Speed** MPH 5 Maximum Climbing Angle/ %(°) 14(8) Safe Climbing Angle<sup>4</sup> Range (Per Charge)<sup>6</sup> Mile 25 **Turning Radius** 50" in **Ground Clearance** 3″ in Handle Bar Type T-Bar Model # 888B-4 **Standard Seat** WxDxBH SV191821B

**Note:** <sup>1</sup> On level Road. <sup>2</sup> Without basket and seating. <sup>3</sup> With basket, seat & batteries. <sup>4</sup> Base on user weight of 200LB. <sup>5</sup> Based on user weight of 165LB. <sup>6</sup> After the batteries and mechanical moving parts fully break in.

