

## THE MOST COMFORTABLE RECLINER IN THE WORLD!







The Voya Collection by **VivaLift!**® brings you superior relaxation and comfort. As you enjoy a refreshing drink, press a button on the remote to adjust the power lumbar or power headrest for complete shoulder support. The Voya recliner's remote also offers the convenience of USB charging. Experience the appeal of the Voya Collection.

**Price Mobility**® (US) 800-800-8586 | (CAN) 888-570-1113 pridemobility.com | prideliftchairs.com | ❸ ② ⑤





COLLECTION	VOYA v.2
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19.5″
SEAT DEPTH	21.5″
SEAT WIDTH	22″
TOP OF BACK TO SEAT	30″
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	16"
OVERALL WIDTH	35″
CHAIR WEIGHT	153 lbs.
WARRANTY	7-Yr. Pro-rated



## **FUNCTIONS AND FEATURES:**

POWER HEADREST



**POWER LUMBAR** 



USB REMOTE



LITHIUM BATTERY BACKUP



FOOTREST EXTENSION



POWER BACKREST

## **POSITION GUIDE:**



RELAXED



RECLINED



**UPWARD LIFT** 



INFINITE LAY FLAT

## STANDARD FABRIC: PICASSO





STEEL

**MIDNIGHT** 

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.