

THE MOST COMFORTABLE RECLINER IN THE WORLD!



METRO



Contemporary comfort has never been more stylish! The Metro Collection by **VivaLift!**® Power Recliners offers maximum relaxation with infinite, lay-flat positioning. Soft, dark colors with Saville fabric offers a unique and cozy appeal, making the Metro recliner a perfect accent piece for the bright areas in your home. A USB port in the remote offers convenient access for charging smart devices. It's time to experience the richness and luxury of the VivaLift!® brand.

Price Mobility® (US) 800-800-8586 | (CAN) 888-570-1113 pridemobility.com | prideliftchairs.com | ❸ ② ⑤





COLLECTION	METRO v.1
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19"
SEAT DEPTH	21″
SEAT WIDTH	21″
TOP OF BACK TO SEAT	26.5″
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	41"
DISTANCE FROM WALL	14"
OVERALL WIDTH	33″
CHAIR WEIGHT	142 lbs.
WARRANTY	7-Yr. Pro-rated



FUNCTIONS AND FEATURES:

POWER HEADREST



POWER LUMBAR



USB REMOTE



LITHIUM BATTERY BACKUP



FOOTREST EXTENSION



POWER BACKREST

POSITION GUIDE:



RELAXED



RECLINED



UPWARD LIFT



INFINITE LAY FLAT

STANDARD FABRIC: SAVILLE QUICK SHIP





GREY

BROWN

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.