

Pride® Lifts

Mobility Made Easy

Outlander LP

The **Outlander LP** Exterior Lift Series is comprised of low-profile lifts designed to take your 10" or 14" mid-wheel drive power chairs with you. Low-profile lifts provide increased visibility while the lift is in use or in the stowed position. The design also provides a shock-absorbing "cradle effect" which protects your power chair while the vehicle is in motion. This is a great exterior lift option for sedans, SUVs, mini-vans, full-size vans and pickup trucks. Simple, hitch-based installation requires little to no modification of your vehicle and allows for easy removal.

BASE WEIGHT ¹	76 lbs. for 10" wheels; 75 lbs. for 14" wheels
MAXIMUM SCOOTER/PC OVERALL WEIGHT ²	350 lbs.
MAXIMUM PRODUCT WIDTH	27"
OVERALL PLATFORM LENGTH	19"
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	20"
MOTOR	Sealed 12 Volt DC
WARRANTY	3-Year Transferable
HITCH CLASS ³	II, III, IV
OPTIONS	• Battery pack w/charger • Swing-away adaptor
VEHICLE PLATFORM	Full-size van, Minivan, Truck, SUV and Crossover
LIFTS CURRENT PRODUCTS ^{1,2}	Jazzy® 1113 ATS Jazzy Air™, Jazzy Select® Jazzy Select® 6, Quantum® 610 Q6 Edge®, J6, Jazzy® 600 ES, Jazzy® 614 HD And equivalent competitor products

1. Vehicle dependent.
2. Lifting weight does not include accessories, upgraded seating and/or hardware.
3. Hitch weight carrying capacity varies by vehicle. Call for details.
4. Does not include any accessories or hardware.

To find a compatible lift for your vehicle and mobility product, use our new

LIFT FINDER

pridemobility.com/lift-finder

FEATURES AND BENEFITS

- 350 lbs. lifting capacity
- Power chair hold-down bar
- Compatible with Class II, III or IV hitches
- Compatible with 10-inch and 14-inch mid-wheel drive power chairs
- Low-profile platform provides better visibility
- Standard license plate holder
- Manual override
- Standard hardware kit



MADE IN
USA

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

