



TranSport Aluminum Transport Chair





FIGURE A

RETAIL BOX

- Weighs 17.5 lbs without footrests
- Aluminum frame is lightweight and strong
- Frame magnets keep footrests secure against frame to prevent swinging open when the chair is carried
- Padded flip-back arms are removable
- Upholstery is comfortable, lightweight and easy-to-clean
- Two-way back release allows to push or pull lever
- The transport chair conveniently folds for travel and storage (Figure A)
- Back folds forward
- Comes with tool-free height adjustable swing-away footrests
- Aluminum rear wheel locks
- Weight Capacity: 250 lbs

PRODUCTS	Description
TS19	Tran Sport Chair FlipBack-RemArms, Blue Upholstery, Blue Frame, 1/cs
CDECIFICATION C	
SPECIFICATIONS	Description
Overall Width	24.5" (W) (open) 10" (W) (closed)
Overall Height	36"
Seat-To-Floor Height	18"
Back Height	16"
Seat Dimensions	19" (W) x 16" (D)
Overall Length w/ Riggings	38"
Weight (without Footrests with Footrests)	17.5 lbs. 20 lbs.
Carton Shipping Weight	19.5 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime
PARTS	Description
STDS871	Handgrip, 1/ea
STDS1A2308	Back Upholstery, 1/ea
STDS1A2311	Arm Release, 1/ea
STDS001013	Folding Cane Hinge, 1/ea
TRB-1	Back Release Lever, 1/ea
STDS1A2307	Seat Upholstery, 1/ea
STDS2A184R	Brake Assm. (Right), 1/ea
STDS2A184L	Brake Assm. (Left), 1/ea
STDS1011	Caster Bearing, 1/ea
STDS2A1811	Rear Axle Bolt, 1/pr
STDS2A1801	Rear Hub Cap, 1/pr



PARTS	Description
STDS2A1800	Wheel, 1/ea
STDS1A2313	Magnet, 1/pr
STDS2A1802	Tipping Anchor, 1/ea
STDS2A1806	Front Axle Bolt, 1/pr
STDS1A2300	Wheel & Fork, 1/ea
STDS2A1810	Front Wheel, 1/ea
STDSA2309R	Footrest- pair TS19SF (Rt), 1/ea
STDSA2309L	Footrest- pair TS19SF (Lt), 1/ea
STDS1A2316	Adj. Knob, 1/pr
STDS2A1803	Cross Frame Assist - Blue, 1/ea
STDS2A18X1	Cross Frame: Blue, 1/ea
STDS1060	Fork Stem Bearing, 1/ea
STDS1A2302R	Side Frame (Rt), 1/ea
STDS1A2302L	Side Frame (Lt), 1/ea
STDS2A1807	Seat Rail Pad, 1/ea
STDS1A2305R	Back Cane & Hinge (Rt), 1/ea
STDS1A2305L	Back Cane & Hinge (Lt), 1/ea
stds1a2318I	Brake Assembly, 1/ea
STDS1A2318R	Brake Assembly, 1/ea
STDS1A2309L	LeftFootrestAssembly, 1/ea
STDS1A2309R	RightFootrestAssembly, 1/ea
TS19SF	FootRest, 1/pr
STDS1069-4	Hanger Pin Set (4), 4/pkg