

Deluxe Fly-Weight Aluminum Transport Chair



FIGURE A



FIGURE B



INSET C

- Aluminum frame weighs only 19 lbs
- Tan Plaid (shown) or Burgundy upholstery (Figure A)
- Composite wheels are lightweight and maintenance free
- Deluxe all-aluminum rear wheel locks
- Comes with carry pocket on backrest
- Both 17" and 19" come with tool-free, swing-away footrests
- Comes with seat belt for added safety
- Back folds down for easy storage and transport (Figure B - C)
- Weight Capacity: 300 lbs

PRODUCTS

PRODUCTS	Description
FW17GR	Green Frame and Black Upholstery, FlyWeight Transport Chair 17Green, 1/cs
FW17DB	Black Frame and Tan Plaid Upholstery, Fly Weight Chr.17TanPlaid, 1/cs
FW17BG	Green Frame and Burgundy Upholstery, Flyweight TransportChr 17Burg/Grn, 1/cs
FW19BL	Blue Frame and Black Upholstery, FlyWeightChair19 AlumBlue, 1/cs
FW19GR	Green Frame and Black Upholstery, Fly Weight Chair, 1/cs
FW19DB	Black Frame and Tan Plaid Upholstery, Fly Weight Chr.19TanPlaid, 1/cs
FW19BG	Green Frame and Burgundy Upholstery, FlyWeight Transport Chr 19Burg/Gr, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Overall Width (FW17)	20" (open) 9" (closed)
Overall Width (FW19)	22" (open) 9" (closed)
Product Overall Height (FW17 FW19)	36" 37.25"
Seat-To-Floor Height	19"
Back Height	18"
Seat Dimensions (FW17)	17" (W) x 16" (D)
Seat Dimensions (FW19)	19" (W) x 16" (D)
Overall Length w/ Riggings	33"
Weight	19 lbs.
Carton Shipping Weight	24 lbs.

ACCESSORIES

ACCESSORIES	Description
STDS1040S	Universal Cup Holder, 1/ea

PARTS

PARTS	Description
TR-17B	back upholstery for s/n with A, 1/ea
STDS1M2017	for s/n with M, 1/ea
STDS001021B	blue backrest cane, 1/ea
TRHG	Handgrip, 1/bx

PARTS	Description
TRB-1	Back Release Lever, 1/ea
SDTS001011	folding cane hinge, 1/ea
SDTS2005	full arm pad nylon, 1/ea
SDTS2004	desk arm pad nylon, 1/ea
STDS002001	arm pad bolts, 1/ea
TR-17S	seat upholstery for s/n with A, 1/ea
FWBC	Clamp for brake, 1/ea
TR-AXLE	8" Rear Wheel Axle, 1/ea
STDS002050	8" Rear Wheel Bearing, 1/ea
STDS002051	rear axle spacer, 1/ea
TRW-1	rear wheel, 1/ea
STDS001044	caster fork, 1/ea
STDS1002-1	Caster axle bolt, 1/ea
STDS1011	Caster Bearing, 1/ea
STDS1066	Dust Cap, 1/ea
STDS4020	tipping sleeve, 1/ea
STDS1002	Standard 8" Caster, 1/ea
STDS1060	Fork Stem Bearing, 1/ea
STDS001045	Fork Stem Nut, 1/ea
FW-1001	front seat rail guide, 1/ea
STDS1070	Seat Rail Guide, 1/ea
STDELR-1L	release lever- Left, 1/ea
STDELR-1R	release lever- Right, 1/ea
STDS001051	Footrest adjustment bolt, 1/ea
STDS1026	Footrest Bumper, 1/ea
FWFP-R	Footplate- right, 1/ea
STDS001061	cross brace bolt, 1/ea
FWSFB	Footrest set: blue, 1/cv
FWSFS	Footrest set: silver, 1/cv
FWSFG	Footrest set: green, 1/bx
FWSFBK	Footrest set: black, 1/ea
STDS2385BG	Back Uphol.17"x18", 1/ea
STDS2385DB	Back Uphol, 1/ea
STDS001021R	Backrest Cane/FW, 1/ea
STDS2A1300	Backrest Cane/FWDB, 1/ea
STDS001011	Folding Cane Hinge FW Series, 1/ea
FWDB-AP	Arm Pad, 1/ea
STDS2A1401	Arm Pad Full Lngth/FWBG, 1/ea
STDSFLNB	Arm Pad- Full, 1/ea
STDS2384DB	Seat Uphol 17"x16", 1/ea
STDS2A1427	Seat Uphol./FWBG, 1/ea
STDS2A4335	Brake Tip, 1/ea
STDS4017L	Brake Assmbly, 1/ea
STDS4017R	Brake Assmbly, 1/ea
STDS1069-4	Hanger Pin Set (4), 4/pkg
STDS5D1070	Seat Rail Guide, 1/ea
FWFP-L	Footplate- left, 1/ea
STDS5D5504	Back Upholstery Screws (8), 1/pkg
STDS2349BLK	Back Upholstery 17", 1/bx
STDS2351BLK	back upholstery for s/n with A, 1/bx
STDS4S3601	Back Cane, 1/ea
STDS4S1503	Folding Hinge TR39, 1/ea