



Duet Rollator/Transport Chair, 8" Casters Padded Seat, Loop Locks



FIGURE A

FIGURE B

FIGURE C

FIGURE D



FIGURE F



795BU



LIFESTYLE

- Combines the features of a rollator and transport chair in one unit An individual can ambulate independently or be pushed safely by a caregiver
- To use as a transport chair the Duet has locking flip-up armrests (Figure A - Figure C) and a flip-down footrest (Figure D) to aid the user A safety lock release is used to store the arms and a hook and loop fastener strap is used to store the footrest while using the Duet for ambulation
- The tool-free two-position contoured padded backrest is reversible depending on whether you are using the Duet as a rollator or in transport chair mode
- Comes with padded seat with carry pouch for added privacy and security of personal items Inside the pouch is a handle to aid with the folding
- The folding frame has a lock-out feature to ensure the unit stays open while in use and ensures the unit will stay folded during transport
- 8" casters are ideal for indoor and outdoor use
- Handles are adjustable in height, with self-threading knob (Figure E)

| PRODUCTS | Description |
|----------|----------------|
| 795B | Blue, 1/cs |
| 795BK | Black, 1/cs |
| 795BU | Burgundy, 1/cs |

| SPECIFICATIONS | Description |
|------------------------|-------------------------------|
| Overall Dimensions | 25.5" (W) x 37" (H) x 26" (L) |
| Seat To-Floor Height | 21" |
| Handle Height | 31.5" - 37" |
| Back Height | 35" |
| Seat Dimensions | 13.5" (W) x 13" (D) |
| Weight | 20 lbs. |
| Carton Shipping Weight | 25 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |

| PARTS | Description |
|--------------|--------------------------------------|
| 9501A79511 | Handle Bar 1A & 2Y ser# Left, 1/ea |
| 9501A79512 | Handle Bar 1A & 2Y ser # Right, 1/ea |
| 750HBAR-2M | Handlebar 2M ser#, 1/ea |
| 950340079901 | Handgrip- Rt., 1/ea |
| 950340089901 | Handgrip- Lt., 1/ea |
| 9501A79501 | Backrest, 1/ea |
| 750BRAKE | Right and left, 1/ea |
| 750BRAKE-2M | Brake 2m ser #, 1/ea |

| PARTS | Description |
|---------------|---|
| 750BCL | Brake Assy. 1A & 2Y ser # left, 1/ea |
| 750BCR | Brake assy. 1A & 2Y ser # right, 1/ea |
| 750BCL-2M | Brake Assy. 2M ser # left, 1/ea |
| 10208KN | Height Adjustment Knob, 1/ea |
| 9501A7951L | Elbow Rest Assembly 1A & 2Y ser # left, 1/ea |
| 9501A7951R | Elbow Rest Assembly 1A & 2Y ser # Right, 1/ea |
| 9502M7951L | Elbow Rest Assm. 2M ser # Left, 1/ea |
| 9502M7951R | Elbow Rest Assm. 2M ser # right, 1/ea |
| 9501A79502 | Folding Hinge Assm. Includes rod and both side bars, 1/ea |
| 9501A7502L | Wheel Brake, Lt., 1/ea |
| 9501A7502R | Wheel Brake, Rt., 1/ea |
| 9500001232 | a. Cable Adjuster, 1/ea |
| 9501A7503 | b. Spring Boot, 1/ea |
| 9500004101 | c. Cable Fastener, 1/ea |
| 951021820000 | d. Cable Spring, 1/ea |
| 95017953R | Rear Wheel, 1/ea |
| 9501A75005 | Rear Axle Cap, 1/ea |
| 9500075009 | Rear Axle Nut, 1/ea |
| STDS1011 | Caster Bearing, 1/ea |
| 9501A75007 | Front Wheel, 1/ea |
| 750FORK | Fork 1A @ 2Y ser #, 1/ea |
| 750FORK-2M | Fork 2M ser #, 1/ea |
| 9501A7953R | Footrest Assembly 1A & 2Y ser # Red, 1/ea |
| 9501A7953B | Footrest Assembly 1A & 2Y ser # blue, 1/ea |
| 9501A7953BK | Footrest Assembly 1A & 2Y ser # black, 1/ea |
| 9501A7953R-2M | Footrest Assembly 2M ser # Red, 1/ea |
| 9501A7953B-2M | Footrest Assembly 2M ser # Blue, 1/ea |
| 795SEAT | Seat, 1/ea |
| 9501A79505 | Elbow Pad, 1/ea |
| 9500001232 | Bushing, 1/ea |
| 950179504 | Folding Hinge Bolt, 1/ea |
| 9501A79504 | Tote Bag 1A & 2Y ser #, 1/ea |
| 9501A79504-2M | Fits Underneath Seat, 1/ea |
| 9501A75913 | Seat Bracket 1A & 2Y ser#, 1/ea |
| 9501A79513-2M | Seat Bracket 2M ser# |
| 9501A79510 | Fork Fastener 1A & 2Y ser # with set screw, 1/ea |
| 9501A79510-RP | Fork Fastener 1A & 2Y ser # with roll pin, 1/ea |
| 9501A79510-2M | Fork Fastener 2M ser #, 1/ea |
| 9501A795010 | Set Screw |
| 9500075010 | Front Axle, 1/ea |
| 9501A79506 | Rear Axle Assembly 1A & 2Y ser #, 1/ea |
| 9501A79506-2M | Rear Axle Assembly 2M ser #, 1/ea |
| 9502M75013 | Sear Bracket 2M, 1/ea |
| 9501A75003 | b. Spring Boot, 1/ea |
| 9501A75008 | 3 Hole Rear, 1/ea |
| 950300169901 | Anti-rattler, 1/ea |
| 9500075005 | Folding Hinge Bolt, 1/ea |
| 9501A75013 | Seat Bracket 1A and 2Y, 1/ea |