



Viper Wheelchair
Deluxe High-Strength, Lightweight, Dual-Axle

HCPCS: K0004




FIGURE A

FIGURE B

- Viper weighs 32 lbs (Excluding front rigging)
- Easier to propel and transport
- Frame style eliminates seat guides and allows for custom backs and accessories
- Back height adjustable from 18"– 20"
- Carbon steel black powder-coated frame double coated for an attractive, chip-proof, maintainable finish
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Composite Mag-style wheels with aluminum hand rims are lightweight and maintenance free
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight
- 8" Front casters have three height adjustments
- Nylon, padded armrests provide added patient comfort
- Flip-back arms allow for easy transfer
- Composite steel lined front receiver is angled so the arm will lock into place every time
- The front mounted release lever can be flipped up or down to unlock flip-up arm, depending on the patient's needs (Figure A - B)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Dual axle provides easy transition of seat height to hemi level
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite foot plates with heel loops are durable and lightweight
- Comes standard with flip-up anti-tippers
- Comes with push-to-lock wheel locks
- Available in 12" and 14" seat widths

PRODUCTS	Description
L412DDA-SF	Viper 12FlpBk.DskArm SF, 1/cs
L412DDA-ELR	Viper 12FlpBk.DskArm ELR, 1/cs
L414DDA-SF	Viper 14FlpBk.DskArm SF, 1/cs
L414DDA-ELR	Viper 14FlpBk.DeskArm ELR, 1/cs

SPECIFICATIONS	Description
Width Open (12" Seat 14" Seat)	20" 22.5"
Width Open (16" Seat 18" Seat)	22.5" 24.5"
Width Open with 20" Seat	28"
Width Closed	12"
Height with 12" and 14" Seats	32"
Height with 16", 18" and 20"Seat	36"
Seat-To-Floor Height (12" - 14") (16" - 20")	13.5" - 15.5" 17.5" - 19.5"
Back Height	16" - 18"
Seat Depth with 12" Seat	12"
Seat Depth with 14" Seat	14"

SPECIFICATIONS	Description
Seat Depth with 16", 18" and 20" Seat	16"
Length with Riggings with 12" and 14" Seats	36"
Length with Riggings with 16", 18" and 20" Seat	42"
Product Weight with 12" Seat	30 lbs.
Product Weight with 14" Seat	31 lbs.
Product Weight with 16" Seat	31.5 lbs.
Product Weight with 18" Seat	32 lbs.
Product Weight with 20" Seat	34 lbs.
Carton Shipping Weight	39 lbs. or less
Weight Limit (with 12" and 14" seat)	250 lbs.
Weight Limit (with 16", 18" and 20" seat)	300 lbs.
Warranty	Limited Lifetime
Warranty	Limited Lifetime
Warranty	Limited Lifetime

ACCESSORIES	Description
STDS1001	Standard 24" Wheel, 1/ea
STDS801	Wheel Lock Extension, 1/ea
STDS1002	Standard 8" Caster, 1/ea
STDS807	Anti-Tippers With Wheels, 1/pr
STDS814	Anti-Tippers With Wheels, 1/pr
STDS840	Flat-free Wheel With Handrims, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS1022	Oblique Projection Handrim, 1/ea
STDS4055L	Left (For use with 18" and 20"), 1/ea
STDS4055R	Right (For use with 18" and 20"), 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS1818N	Black Extension, 1/ea
STDS2018N	Seat Rail Extension Kit, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea

PARTS	Description
STDS1001	Standard 24" Wheel, 1/ea
STDS801	Wheel Lock Extension, 1/ea
STDS1002	Standard 8" Caster, 1/ea
STDS826	Amputee Adapter, 1/pr
STDS826C	Amputee Adapter, 1/pr
STDS807	Anti-Tippers With Wheels, 1/pr
STDS814	Anti-Tippers With Wheels, 1/pr
STDS835	1/ea
STDS840	Flat-free Wheel With Handrims, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS1022	Oblique Projection Handrim, 1/ea
STDS4055L	Left (For use with 18" and 20"), 1/ea
STDS4055R	Right (For use with 18" and 20"), 1/ea